
INORGANICS – CHEMICAL INDUSTRY

APPLICATION SHEET

CALCIUM HYPOCHLORITE

Calcium hypochlorite $Ca(OCI)_2$ is an oxidizing substance and is used i.e. in detergents as bleech. It is also used as disinfection material of drinking water containers and swimming pools. As calcium hypochloriteis an oxidizer it likes to react with organic materialswhen in direct contact and heated.



Instrument TG 209 F1 – QMS Aëolos®

Test Conditions

Temperature range RT ... 600°C Heating rate 10 K/min Atmosphere Nitrogen at 40 ml/min Sample mass 18.9 mg Crucible Alumina Sensor TG type Platinel

Results

Calcium hypoclorite was measured with a TG with QMS system connected to the gas outlet. The original "water-free" calcium hypochlorite had hydrolysed during storage. This can clearly be seen from the TG and mass spectrometer results during the 1st, 2nd and 4th TG mass losses where water is detected (84°C, 189°C and 489°C). Pure "water-free" calcium hypochlorite shows only two TG steps where oxygen is released (2nd and 3rd TG steps). The remaining product would be CaCl₂ at 400°C. Because of the water content, also small amounts of HCl were found (74°C and 197°C, peak temperatures).

