APPLICATION SHEET



## POLYMERS – POLYMER MANUFACTURING

# POLYCARBONATE (PC)

Polycarbonates are a particular group of thermoplastics. They are easily worked, molded, and thermoformed. They are called polycarbonates because these polymers have functional groups linked together by carbonate groups (-O-CO-O-) in a long molecular chain. The most common type of polycarbonate plastic is one made from Bisphenol A, in which groups from Bisphenol A are linked together by carbonate groups in a polymer chain. Polycarbonate is becoming more common in house wares as well as laboratories and industry. It is often used to create protective features, for example, in banks as well as vandal-proof windows and lighting lenses for many buildings. Other products made from polycarbonate include sunglass/eyeglass lenses, compact discs, DVDs, and automotive headlamp lenses.



### Instrument TG 209 F1 Iris<sup>®</sup> – QMS 403 Aeolos<sup>®</sup>

### **Test Conditions**

Temperature range Heating/cooling rates Atmosphere Sample mass Crucible Sensor RT ... 700°C 20 K/min Nitrogen at 40 ml/min 7.3 mg Alumina Platinel

#### Results

During the decomposition of PC under inert conditions, the polymer chains, of course, decompose and fragments of Bisphenol can be detected indicated i.e. by mass numbers 77 and 78, representing benzene. But also higher mass numbers occur, for example, 120 and 122 amu.



NETZSCH