
APPLICATION SHEET





Motor oil is used as a lubricant in various kinds of internal combustion engines in automobiles and other vehicles, boats, lawn mowers, trains, airplanes, etc. In engines, there are contacting parts which move against each other at high speeds, often for prolonged periods of time. Such rubbing motion causes friction, absorbing otherwise useful power produced by the motor and converting the energy to useless heat. Lubricating oil makes a film

between surfaces of parts moving against each other so as to minimize direct contact between them decreasing friction, wear, and production of excessive heat. Also motor oil carries away heat from moving parts. During usage, engine oil is contaminated with particles. Thermal treatment is also a reason for aging of the oil. In order to avoid engine failure, the oil must be changed when loosing its performance.



Instrument

TG 209 **F3** Tarsus®

Test Conditions

Temperature range Heating rate

Atmosphere Nit Sample mass 10.

Crucible Atmosphere RT ... 1000°C 20 K/min

Nitrogen at 20 ml/min to 700°C

10.14 and 11.56 mg Al₂O₃

 Ai_2O_3 air to 1000°C

Results

The soot content of engine oil is a good indicator of the maintenance interval. The main decomposition step between 200 and 450°C is related to the pyrolysis of the oil. At 700°C, the purge gas is changed from nitrogen to air. If carbon is present, it reacts with oxygen and forms CO₂. The carbon dioxide release can be determined as a mass-loss step. The figure shows no carbon content for the new oil whereas the used oil contains 6.5% soot. For soot determination, optical methods were formely used. However, other particles from metal abrasion for instance often disturb these analyses or increase the uncertainty.

