

Temperature Dependance of Polyetherimide (PEI)

This figure shows the course of the complex elasticity modulus and the visco-elastic damping $tan\delta$ in dependence of the temperature. The marked modulus decline for about 4 decades during the glass transition, starting at 215°C, is

particularly significant. Here, the material begins to flow. From the low-temperature region (-150°C) up to temperatures above the glass transition, damping increases by 3 orders of magnitude. Material changes from solid to liquid are detectable by this technique.



Testing of a Polymer Blend

continuous decrease in the complex modulus of elasticity with temperature. The line shape of the tan δ dependence on the temperature shows phase transitions at -90°C, 0°C and 100°C.







1 1 NETZSCH-Gerätebau GmbH Wittelsbacherstraße 42 · 95100 Selb · Germany Phone: +49 9287/881-0 · Fax: +49 9287/881505 at@netzsch.com · www.netzsch.com