APPLICATION SHEET

Metals – DMA EPLEXOR®

NiTi Shape Memory Alloys – New Materials Investigated with the Powerful DMA Method

This plot shows the complex modulus E^* and loss factor tan δ of an TiNi wire (shape memory alloy) as a function of temperature. With increasing temperature (heating rate:

3 K/min, test frequency: 10 Hz, tensile mode), the shape memory alloy undergoes the transformation from martensite to austenite at a material specific transition temperature of 100°C. At the transition temperature, a sudden increase in tensile modulus E^* as well as a decrease in intrinsic visco-elastic damping tan δ can be observed.

Temperature Dependence of Magnesium

The figure shows the complex modulus E^* and loss factor tan δ of a magnesium bar (thickness: approx. 2 mm, width: approx. 4 mm) as a function of temperature. With increasing temperature (heating rate: 2 K/min, test frequency:

10 Hz, modified single cantilever), the modulus decreases softly about 20%. The intrinsic damping shows a significant increase. A strong increase for tan δ of about four orders of magnitude has to be noticed. This run shows the very high resolution of the tan δ measurement due to the very high resolution and accuracy of the EPLEXOR[®].





1 1 NETZSCH-Gerätebau GmbH Wittelsbacherstraße 42 · 95100 Selb · Germany Phone: +49 9287/881-0 · Fax: +49 9287/881505 at@netzsch.com · www.netzsch.com