APPLICATION SHEET




## POLYMERS – POLYMER MANUFACTURING

# POLYPROPYLENE

Polypropylene (PP) is a thermoplastic polymer, used in a wide variety of applications, including food packaging, textiles, laboratory equipment, automotive components,

and polymer cash cards. As an additional polymer made from the monomer propylene, it is unusually resistant to many chemical solvents, bases and acids.



#### Instrument

DSC 200 F3 Maia® / TG 209 F3 Tarsus®

### **Test Conditions**

```
Temperature range
Heating rate
```

Atmosphere Sample mass (TG) Sample mass (DSC)  $25^{\circ}$ C ...  $200^{\circ}$ C / N<sub>2</sub> 200°C ...  $300^{\circ}$ C / O<sub>2</sub> 5 K/min Nitrogen at 20 ml/min Air at 20 ml/min 11.93 mg 12.95 mg

#### Results

The measurements were carried out with the TG and DSC systems. The DSC shows the melting peak at 166°C and the onset of degradation at 22.9°C. In the TG run, the c-DTA signal can be calculated during the measurement. Therefore, the O.I.T. (Oxidation Induction Temperature) values can be measured with the DSC and also with the thermobalance. The O.I.T. values are given by the onset of the exothermic degradation peak of the DSC or c-DTA curve. They are comparable between the DSC test (222.9°C) and the thermogravimetric test (222.8°C).





NGB