




## Investigation of a PP-CNT Compound by Means of *Identify*

The new *Identify* software which is an extension of the NETZSCH *Proteus®* software ...

- ... is a unique DSC curve recognition and interpretation system providing results with a single click.
- ... is useful for material identification and quality control.
- ... is both easy to use and sophisticated
- ... includes a database with NETZSCH libraries for polymers as a basis as well as libraries that can be created by the user.
- ... manages measurements, literature data and classes, incorporating the user's knowledge.



DSC measurement (2<sup>nd</sup> heating) on a PP+NT compound investigated by means of *Identify* 

A PP+CNT (polypropylene + carbon nanotubes) compound was investigated by means of a DSC measurement (2<sup>nd</sup> heating) using *Identify*. The highest similarity of 48.6% between the measured DSC curve (white) and database entries can be found for the literature data for PP\_isotactic included in the NETZSCH "Thermal Properties of Polymers" poster.

Although the DSC curve was measured at a relatively high heating rate of 40 K/min, PP could be identified correctly. A database measurement on PP measured at 10 K/min (also  $2^{nd}$  heating) is displayed in pink for comparison.

