Technical Specifications



STA 2500 Regulus	
Top-loading, vertical & differential balance	
RT to 1100°C/1600°C (2 easily exchangeable furnaces)	
0.001 to 100 K/min (depending on the furnace)	
0.3 K	
Single hoist	
Differential balance system	
± 250 mg, max. sample load 1 g	
0.03 μg	
Type S	
10^{-4} mbar (10^{-2} Pa), depending on the evacuation system	
Inert, oxidizing, vacuum	
Built-in mass flow controllers	
Optionally available for 20 crucibles	